

## How to pair your hearing aids with your Android™ device

1. Go to **Google Play**. In the **Search** field on top, enter the name of the app and download the app.
2. Once downloaded, open the app and tap **Get Started**.
3. Follow the instructions on the screen to connect. You are now ready to go.

**NOTE:** If you have previously paired any hearing aids, you need to unpair them in your smartphone's Bluetooth menu before you pair them again.

## If you need to re-pair your iOS device

1. Start by removing the old pairing by going to the **Settings** menu.
2. Tap **General**, scroll down and select **Accessibility**.
3. Select **MF Hearing Devices** and tap on **Forget this device**.
4. Turn your iOS device off and on again and re-pair the hearing aids to the device. If you need more information, contact your hearing care professional.

Interton GN

### Worldwide Headquarters

Interton A/S  
Lautrupbjerg 7  
DK-2750 Ballerup  
Denmark  
Tel.: +45 4575 1111  
interton.com

CVR no. 55082715

### United Kingdom

GN Hearing Ltd.  
Kirtlington Business Centre  
Portway  
Kirtlington  
Oxon OX5 3JA  
Tel.: +44 1869 352 800  
interton.com

### Australia

GN Hearing Australia Pty Ltd  
Gate C, 19-25 Khartoum Road  
Macquarie Technology Park  
Macquarie Park NSW 2113  
Tel.: (free) 1800 658 955  
interton.com

### Singapore

GN Hearing Pte. Ltd.  
2 Kallang Avenue  
#07-19 CT HUB  
Singapore - 339407  
Tel.: +65 6320 9388  
interton.com

### New Zealand


GN Hearing New Zealand Limited  
Ground Floor, North Entrance  
4 Fred Thomas Drive  
Takapuna, Auckland, 0622  
Tel.: (free) 0800 900 126  
interton.com

© 2019 GN Hearing GmbH. All rights reserved. Interton is a trademark of GN Hearing GmbH. Apple, the Apple logo, iPhone, iPad, iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC. The Bluetooth word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.

Interton GN

## Get connected

The Interton Sound app



M201465GB-19.11-Rev.B

GN Making Life Sound Better  
FOR 150 YEARS

Made for  
iPhone | iPad | iPod

## Introduction


Thank you for choosing our products. We have created this app to make it easy and pleasant for you to use and control your hearing aids from your smartphone or mobile device.

For information on which smartphones and mobile devices are compatible with the app, visit: [intertron.com/compatibility](http://intertron.com/compatibility)

Before you can start using the app, you must pair your smartphone or mobile device with your hearing aids. This is easily done - just follow the steps below.


## How to pair your hearing aids with your iOS device

1. Turn your hearing aids off.
2. Turn on **Bluetooth** by going to the **Settings** menu on your iOS device. The slider turns green when you activate Bluetooth.



3


3. Go back to the main **Settings** menu and then, tap on **Accessibility**.
4. Scroll down to **Hearing Devices**.



4

5. Turn your hearing aids on. Now, your iOS device will be able to detect your hearing aids.

**NOTE:** From iOS 13 and newer.



6

6. Tap on **Pair**.
7. Go to the **App Store**, and in the **Search** field on top enter the name of the app.
8. Download the app, open it and tap **Get Started**.
9. Follow the instructions on the screen to connect to the app. You are now ready to go.

