

LIVE LIFE FROM BOTH SIDES – DIGITAL WIRELESS CROS

WIRE & ADAPTER
FREE




) RELAY (

USER INFORMATION


LIVING IN STEREO

Having two functioning ears is extremely important in everyday life. When one ear suffers from a total loss of hearing, many day to day situations become difficult to manage. Understanding each member of a group around the table, walking and talking side by side or even crossing a road become an utterly exhausting task.

CROS Hearing Aids (respectively BiCROS for cases where even the better ear needs additional help) have been known to be the best solution when the hearing in one ear is much worse than the other. CROS stands for "Contralateral Routing Of Sounds". Contralateral means "to the other side". Simply said a CROS hearing aid receives the sound arriving on one side of your head and sends it into your opposite ear.



- **Poorer ear:**
No hearing ability
- **Better ear:**
Normal hearing



- **Poorer ear:**
No hearing ability
- ◐ **Better ear:**
Hearing loss



FREEDOM, FREE OF INTERFERENCES

Traditional CROS devices have been around for a while. Until now CROS solutions suffered from some inconvenient flaws. Disturbing interference from radio or cell phones, oversized and uncomfortable adapters or unattractive wires led to less than satisfactory results.

Dedicated to solve these issues, Interton has developed RELAY, the first digital wireless “all-in-one-CROS-System” in the world. RELAY consists of one transmitter and one receiver unit to fit to each ear.

By using digital wireless transmission RELAY is the one wireless CROS system not requiring an additional adapter for the sound transfer. So you can really have it all, the freedom of a wireless system and interference-free digital sound quality.


RELAY - making CROS hearing better and more comfortable than ever.

REACH THE OTHER SIDE!



RELAY – DIGITAL WIRELESS CROS SYSTEM

- Advanced digital hearing technology
- Better understanding of speech arriving at the bad ear
- Available in "Open-Fit" (for extra comfort), Standard ear-hook (for custom fit) and In-The-Ear versions
- "All-In-One" System, no adapters or cables needed
- Interference-free digital wireless transmission
- For CROS, BiCROS and monaural fittings
- Improved awareness of surrounding sound
- Better hearing even when the good ear needs assistance
- The smallest wireless CROS system (Height: 40 mm, Width: 8,6 mm)




Solutions available:



Thin tube/Open-Fit



Standard Ear-Hook



In-The-Ear

Colours available:



beige



brown



grey



Steelgrey



Customer contact

Interton Electronic Hörgeräte GmbH | Am Dännekamp 15 | D-51469 Bergisch Gladbach | Germany
Tel.: +49 2202 95 26 31

Worldwide headquarters

Interton A/S | Lautrupbjerg 9 | P.O. Box 130 | DK-2750 Ballerup | Denmark
Tel.: +45 45 75 11 11 | Fax: +45 45 75 11 19 | www.interton.com