











## LISTENING SITUATIONS:

These are typical everyday scenarios where you might find it difficult to e.g. follow conversations. To the right you find an indication how **AVIO1** suits the respective listening situations.

| <b>avio<sup>1</sup></b>   |           |
|---|-----------|
|  Telecoil Situation<br>· On the phone / at church   | ● ● ● ● ● |
|  Speech in Quiet<br>· 1:1 Conversations at home   | ● ● ● ● ● |
|  Situations without Speech<br>· A car passing by<br>· Walking noises                      | ● ● ● ● ● |
|  Passive Music<br>· TV / radio / home stereo  | ● ● ● ● ● |
|  Live Music<br>· Theatre / opera / concerts  | ● ● ● ● ○ |
|  Loud Environments<br>· Party / Disco<br>· At extreme working places                    | ● ● ● ● ○ |
|  Office Situation<br>· Direct conversation<br>· Printer / computer noise                | ● ● ● ● ○ |
|  Speech in Constant Background Noise<br>· While driving a car                           | ● ● ● ● ○ |
|  Speech in Noisy Situations<br>· Group conversations<br>· Talking outdoors / in traffic | ● ● ● ○ ○ |
|  Always changing Situations<br>· From group meetings to 1:1 conversations               | ● ● ● ○ ○ |

# ALL IT TAKES TO SMILE



# avio<sup>1</sup>

**BASIC CLASS  
HEARING SYSTEM**

#### Customer contact

Interton Electronic Hörgeräte GmbH | Am Dännekamp 15  
 D-51469 Bergisch Gladbach | Germany  
 Tel.: +49 22 02 95 26 31

#### Worldwide headquarters

Interton A/S | Lautrupbjerg 9 | P.O. Box 130  
 DK-2750 Ballerup | Denmark  
 Tel.: +45 45751111 | Fax: +45 45751119 | www.interton.com

M100936-GB-08.12 Rev.A

# audio<sup>1</sup>

BASIC CLASS  
HEARING SYSTEM

## YOUR HEARING AT HEART

Because you want to get the most out of every moment, we put our heart and soul into creating **audio<sup>1</sup>**. Built with your hearing in mind and suitable for most day-to-day situations, **audio<sup>1</sup>** provides basic digital hearing systems at a very attractive price.

Equipped with a set of useful basic functions, **audio<sup>1</sup>** is designed to make your life easier. Numerous enjoyable pursuits are twice as nice when you can communicate with ease.

**audio<sup>1</sup>**'s Adaptive Feedback Cancellation function plays a vital part in this by preventing any unpleasant squealing in your ear, thus helping you to focus on the sounds that really matter.

## YOUR **audio<sup>1</sup>** BENEFITS:

- Trouble-free telephone functionality (Telecoil)
- Up to three personal programs ensure good hearing in various listening situations
- Front focus for better speech recognition
- Background noise reduction decreases distracting noise, which makes relevant sounds clearer
- Feedback Manager for whistle-free hearing experience

Now isn't that something to smile about?

## MODELS

The variety of models available delivers further proof that **audio<sup>1</sup>** is all about ensuring happy hearing.



### Behind the Ear (BTE)

BTE hearing aids rest comfortably behind the ear. Sound is directed from the hearing aid through tubing into the ear canal.



### Open Fit (BTE)

Open Fit hearing systems direct sound from the hearing system through a sleek and discreet thin tube open fit solution. Open fit eliminates the 'plugged up' feeling resulting in a more comfortable listening experience.



### Receiver in the Ear (RIE)

RIE hearing systems are a combination of BTE and ITE models. Sounds are collected by the unit that sits behind the ear. A receiver sits smoothly in the ear canal and transmits the amplified sounds directly to the ear, providing exceptional sound quality.



### In the Ear (ITE)

ITE systems are custom moulded to fit in the ear, providing excellent sound quality and fit. Sizes range from the virtually invisible Completely-In-The-Canal model (CIC) to the Full Shell (FS), which sits in the canal as well as the outer ear.

## **audio<sup>1</sup>** ON EAR

Can you spot the hearing aid? The Behind-the-Ear model is very small.



## COLOURS

**audio<sup>1</sup>** is available in a selection of neutral colours.



Black



Marble Grey



Anthracite



Dark Brown



Beige



Medium Blonde



Caramel

### Additional colours for RIE:



Silver



Glossy White



Glossy Black