








# Cable and wireless fitting guide

Interton Fitting	 <b>20310800</b> (NOAHLINK WIRELESS)	 <b>9022 907 69019</b> (HI-PRO CABLE/RIGHT/CS44) <b>9022 907 69029</b> (HI-PRO CABLE/LEFT/CS44) <b>3122 150 73781</b> (FLEX CONNECTOR CLIP CS63) <b>19181700</b> (STRIP, FLEX, CS63, PRG)	 <b>15206700</b> (PROGRAMMING CABLE CS63, RIGHT) <b>15206701</b> (PROGRAMMING CABLE CS63, LEFT) <b>19181700</b> (STRIP, FLEX, CS63, PRG)	 <b>9022 907 69019</b> (HI-PRO CABLE/RIGHT/CS44) <b>9022 907 69029</b> (HI-PRO CABLE/LEFT/CS44) <b>19533900</b> (ADAPTER, FLEX, 3 PIN, CS53)	 <b>9022 907 69019</b> (HI-PRO CABLE/RIGHT/CS44) <b>9022 907 69029</b> (HI-PRO CABLE/LEFT/CS44) <b>19902500</b> (PROGRAMMING ADAPTOR 2)
Interton Ready					
RIE WL - RDx62	●				●
RIE WL - RDx61	●	●	●	●	
BTE WL - RDx70, RDx80	●				●
Custom WL	●	●	●	●	
Custom non-WL		●	●	●	

**Please note:**

Products supported in Interton Fitting do not support the Airlink dongle. Interton hearing instruments for both wireless and non-wireless fittings **require a fresh battery** when fitted with Interton Fitting.