





CIC

Model	MV6CIC-W	MV4CIC-W	MV3CIC-W
Device Configurations			
Battery size	10A Zinc-Air		
Receiver power levels	LP, MP & HP		
Control option	Push button		
Sound quality			
WARP compression (WDRC), number of channels	17	12	8
Comfort			
Adaptive Noise Reduction	●	●	●
Impulse Noise Reduction	●	●	●
Microphone Noise Reduction	●	●	●
Environmental Gain Tuner	●	●	●
Environmental Classifier	●	●	●
Feedback management			
Feedback Manager Plus	●	●	●
Music Mode	●	●	●
Preset Feedback Manager	●	●	●
Onboarding			
Synchronised Acclimatization Manager	●	●	●
Acclimatization Manager	●	●	●
Convenience			
Ear to Ear Communication (Push Button)	●	●	●
Power-on-delay	●	●	●
AutoPhone	●	●	●
Comfort Phone	●	●	●
Direct audio streaming (MFi, Android™*)	●	●	●
TV Streamer 2, Remote Control, Remote Control 2, Phone Clip 2, Micro Mic and Multi Mic	●	●	●
Interton Sound™ app	●	●	●
Remote Firmware Update	●	●	●
Fitting Features			
Interton Fitting™ 1.10 or higher	●	●	●
Number of Programs	4	4	4
Tinnitus Sound Generator	●	●	●
Datalogging	●	●	●
Wireless Fitting with Noahlink Wireless	●	●	●

* Compatible with Android smartphones that support direct Android streaming to hearing aids.

Fitting range




Technical Specifications


		LP		MP		
		IEC 60118-0:1983_AMD1:1994 IEC 60118-0:2015(*) IEC 711 Ear simulator	ANSI S3.22-2014 IEC 60118-0:2015 JIS C 5512: 2015 2cc coupler	IEC 60118-0:1983_AMD1:1994 IEC 60118-0:2015(*) IEC 711 Ear simulator	ANSI S3.22-2014 IEC 60118-0:2015 JIS C 5512: 2015 2cc coupler	
Reference test gain (60 dB SPL input)	1600 Hz/HFA	33	32	40	37	dB
Full-on gain (50 dB SPL input)	Max. 1600 Hz/HFA	49	40	59	50	dB
Maximum output (90 dB SPL input)	Max. 1600 Hz/HFA	124	114	128	118	dB SPL
Total harmonic distortion	500 Hz	0.5	0.4	0.7	0.8	%
	800 Hz	0.5	0.5	1.1	0.9	
	1600 Hz	0.5	0.7	0.8	1.0	
	3200 Hz	-	0.1	-	0.3	
Equivalent input noise, w/o Noise reduction		22	22	25	24	dB SPL
1/3 Octave Equivalent input noise, w/o Noise reduction	1600 Hz	10	10	11	11	dB SPL
Frequency range IEC 60118-0: 2015		100-8170*	100-7230	100-8250*	100-7970	Hz
Current Drain (Quiescent / Operating)		1.12/1.14	1.12/1.22	1.10/1.13	1.10/1.30	mA
Weight of hearing aid		1.62 / 0.06		1.78 / 0.06		gram/oz

* Measured according to IEC 60118-0:2015, with 711-Ear simulator coupler.



IEC 60118-0:1983_AMD1:1994
IEC 711 Ear Simulator



ANSI S3.22-2014
IEC 60118-0:2015
JIS C 5512: 2015
2cc coupler



■ LP
■ MP



Patents pending


All specifications are subject to change without notice

Technical Specifications


		HP		
		IEC 60118-0:1983_AMD1:1994 IEC 60118-0:2015(*) IEC 711 Ear simulator	ANSI S3.22-2014 IEC 60118-0:2015 JIS C 5512: 2015 2cc coupler	
Reference test gain (60 dB SPL input)	1600 Hz/HFA	47	42	dB
Full-on gain (50 dB SPL input)	Max. 1600 Hz/HFA	69	60	dB
Maximum output (90 dB SPL input)	Max. 1600 Hz/HFA	130	120	dB SPL
Total harmonic distortion	500 Hz	0.5	0.4	%
	800 Hz	1.0	0.8	
	1600 Hz	0.8	0.3	
	3200 Hz	-	0.2	
Equivalent input noise, w/o Noise reduction		25	23	dB SPL
1/3 Octave Equivalent input noise, w/o Noise reduction	1600 Hz	11	11	dB SPL
Frequency range IEC 60118-0: 2015		100-7370*	100-6790	Hz
Current Drain (Quiescent / Operating)		1.17/1.20	1.17/1.24	mA
Weight of hearing aid		1.82 / 0.06		gram/oz

* Measured according to IEC 60118-0:2015, with 711-Ear simulator coupler.



IEC 60118-0:1983_AMD1:1994
IEC 711 Ear Simulator



ANSI S3.22-2014
IEC 60118-0:2015
JIS C 5512: 2015
2cc coupler



■ HP



Patents pending

All specifications are subject to change without notice