Interton. Part of something bigger.

Interton was founded in Germany in 1962 and is today part of GN, a company dedicated to making life sound better through intelligent audio solutions that let you hear more, do more and be more than you ever thought possible.

We are the only company in the world with medical, professional and consumer grade audio technology under one roof. This gives us unprecedented knowledge of hearing and the human ear.

Are you ready to hear better and smile more?

www.interton.com

Interton | Lautrupbjerg 7 | DK-2750 Ballerup Tel: +45 4575 1111 | www.interton.com

© 2021 GN Hearing A/S. All rights reserved. Interton is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC.

d102842GB-21.03-Rev.A



If you like getting more from life, you'll like Interton Move





Life is for living

If you're the kind of person who likes getting more out of life, then it might be time to grab yourself a pair of new Interton Move hearing aids.

At Interton, we've done the thinking about what makes a great hearing aid for you - using GN technology, a world leading heraing aid manufacturer:

- Want to keep track of the conversation, even in noisy environments like your favourite coffee shop or on a windy walk? Of course you do, so Interton Move works together with your own ears to support you no matter the situation.
- Like listening to music or taking your phone calls on the go? So do we - seamlessly stream your favourite tunes, podcasts and calls from your smartphone.
- Want a rechargeable device that goes for as long as you do? Interton Move gives you up to 30 hours of power in as little as 3 hours of charge time.
- Tired of fiddling around with tiny switches? We hear you. Interton Move comes with an easy to use app making them easier to operate.

Put simply, if you think life is for living, we think you'll like Interton Move.

Designed for the life you lead

At Interton, we understand you prefer to get out and enjoy everything life has to offer without having to spend time thinking about your hearing loss or old hearing aids.

This is why we have created Inaterton Move – our clever new hearing aid designed specifically for people who like getting more out of life with clear sound, seamless connectivity and a rechargeable battery option that can keep going as long as you do.

Stay part of the conversation

For people who like socialising and getting out and about, feeling excluded from the conversation simply because you struggle to hear in challenging environments is particularly frustrating.

So, we have upgraded Interton Move with our latest and greatest technology - which now features more computing power, clearer sounds and less noise. This means you can keep track of the conversation, even on a windy walk or in a noisy cafe and have better awareness of the sounds around you.

With Interton Move, getting more out of life means never being sidelined.



Full day of freedom

With Interton Move there's no need to let the fun stop, because - with up to 30 hours of battery power, our new rechargeable solution can keep going for as long as you do. Best of all, go from 0-100% in just 3 hours of charge time.

The Desktop charger is a small, sturdy unit plugged into a wall outlet and it utilizes wireless, inductive charging.

Desktop charger with Interton Move rechargeable hearing aids

It is easy to use, simply remove your hearing aids from your ear and place them in the charger until they are ready to use again. The visual LED indicators on top of the hearing aids keep you informed of the charge status:

- LED slowly flash green: your hearing aids are charging.
- LED with steady green light: your hearing aids have reached a full charge.

Connect to happiness

If getting more out of life means connecting more to the things you love, then we think you'll like Interton Move.

Interton Move seamlessly streams crystal clear music, podcasts or calls direct from your iOS and selected Android[™] devices. And with the Interton Sound app you have control at your fingertips, so you always have the best sound quality to match wherever your day takes you.

You can also enhance your everyday life with our discreet wireless accessories. You can stream crystal clear sound to your hearing aids from your TV, or from a speaker in a oneto-one setting or large gatherings.









With Interton Sound app, you can:

- Make quick, discreet adjustments in real time.
- Locate misplaced devices with the "Find my hearing aid" feature.
- Get the most out of your hearing aids with tips in the app.



Peace of mind

Getting more out life means being able to rely on the things you trust. Feel good knowing Interton Move is a full family of hearing aids packed with GN technology, a world leading hearing aid manufacturer. We give you peace of mind to live life the way you want, rather than having to think about your hearing loss or old hearing aids.

Interton Move comes in all shapes, sizes and colours, so feel confident there's one to suit every ear and every pocket.

Every Interton Move is designed for comfort and durability. The ergonomic design helps you fit and forget, so it stays comfortable all day, while the protective nano-coating makes it resistant to sweat and dust. Life can be tough, wearing your hearing aid shouldn't be.

Designed for a life to be well lived

Receiver-In-the-Ear (RIE) hearing aids





Mini RIE

Rechargeable RIE

Customs hearing aids





Completely-In-the-Canal (CIC-W)

In-the-Canal (ITC)

Behind-the-Ear (BTE) hearing aids



BTE

Power BTE

Interton Move comes in five different colours to match you skin and hair tone.